

method name	used for which samples	pages
HCD_105min_T20_grad3_TB	QuaNPA MAA bead and protein input optimizations	19-22
HCD_2H5M_T20_HI	DDA method for Hela SILAC benchmark	15-18
HCD_105min_DIA_26iw_TB	DIA method 1 for Hela SILAC benchmark	2-5
HCD_105min_DIA_3xhr_MS1_27iw_MS2_TB	DIA method 2 for Hela SILAC benchmark	6-10
HCD_90min_grad3_DIA_28iw_TB	DIA method for IFNg time course NSP samples	11-14

Instrument Method: HCD_105min_DIA_26iw_TB

Thermo EASY-LC method print for QExactive_HF

Sample pickup:

Volume [µl] : 1.00
Flow [µl / min] : 20.00

Sample loading:

Volume [µl] : 22.00
Flow [µl / min] : (unspecified)
Max. pressure [Bar] : 800.00

Gradient:

Time [mm:ss]	Duration [mm:ss]	Flow [nl/min]	Mixture [%B]
00:00	00:00	300.00	3.00
04:00	04:00	300.00	8.00
06:00	02:00	300.00	10.00
74:00	68:00	300.00	32.00
86:00	12:00	300.00	50.00
87:00	01:00	300.00	100.00
94:00	07:00	300.00	100.00
95:00	01:00	300.00	3.00
105:00	10:00	300.00	3.00

Pre-column equilibration:

Volume [µl] : 12.00
Flow [µl / min] : (unspecified)
Max. pressure [Bar] : 800.00

Analytical column equilibration:

Volume [µl] : 3.50
Flow [µl / min] : (unspecified)
Max. pressure [Bar] : 600.00

Autosampler wash:

Flush volume [µl] : 100.00

*Method of Q Exactive HF***Overall method settings****Global Settings**

Use lock masses	best
Lock mass injection	-
Chrom. peak width (FWHM)	20 s

Time

Method duration	105.00 min
-----------------	------------

Customized Tolerances (+/-)

Lock Masses	-
Inclusion	-
Exclusion	-
Neutral Loss	-
Mass Tags	-
Dynamic Exclusion	-

*Experiments***Full MS - SIM****General**

Runtime	0 to 105 min
Polarity	positive
In-source CID	0.0 eV

Full MS - SIM

Microscans	1
Resolution	60,000
AGC target	3e6
Maximum IT	20 ms
Number of scan ranges	1
Scan range	400 to 1000 m/z
Spectrum data type	Profile

DIA**General**

Runtime	0 to 105 min
Polarity	positive
In-source CID	0.0 eV
Default charge state	2

DIA

Microscans	1
Resolution	30,000
AGC target	1e6
Maximum IT	50 ms
Loop count	26
MSX count	1
MSX isochronous ITs	on
Isolation window	23.3 m/z

Isolation offset	0.0 m/z
Fixed first mass	200.0 m/z
(N)CE / stepped (N)CE	nce: 27
Spectrum data type	Profile

Setup**Tunefiles****General**

Switch Count 0
Base Tunefile C:\Xcalibur\methods\nanosource.mstune

Contact Closure**General**

Used False
Start in Closed True
Switch Count 0

Syringe**General**

Used False
Start in OFF True
Stop at end of run False
Switch Count 0

Pump setup

Syringe type Hamilton
Flow rate 3.000 µL/min
Inner diameter 2.303 mm
Volume 250 µL

Divert Valve A**General**

Used False
Start in 1-2 True
Switch Count 0

Divert Valve B**General**

Used False
Start in 1-2 True
Switch Count 0

Lock Masses

1 entry

Mass [m/z]	Polarity	Start [min]	End [min]	Comment
445.12003	Positive			

Inclusion List

26 entries

Mass [m/z]	Formula [M]	Species	CS [z]	Polarity	Start [min]	End (N)CE [min]	MSX ID	Comment
413.65000				Positive				
435.95000				Positive				
458.25000				Positive				
480.55000				Positive				
502.85000				Positive				
525.15000				Positive				
547.45000				Positive				
569.75000				Positive				
592.05000				Positive				
614.35000				Positive				
636.65000				Positive				
658.95000				Positive				
681.25000				Positive				
703.55000				Positive				
725.85000				Positive				
748.15000				Positive				
770.45000				Positive				
792.75000				Positive				
815.05000				Positive				
837.35000				Positive				
859.65000				Positive				
881.95000				Positive				
904.25000				Positive				
926.55000				Positive				
948.85000				Positive				
971.15000				Positive				

Instrument Method: HCD_105min_DIA_3xhr_MS1_27iw_MS2_TB

Thermo EASY-LC method print for QExactive_HF

Sample pickup:

Volume [µl] : 1.00
Flow [µl / min] : 20.00

Sample loading:

Volume [µl] : 22.00
Flow [µl / min] : (unspecified)
Max. pressure [Bar] : 800.00

Gradient:

Time [mm:ss]	Duration [mm:ss]	Flow [nl/min]	Mixture [%B]
00:00	00:00	300.00	3.00
04:00	04:00	300.00	8.00
06:00	02:00	300.00	10.00
74:00	68:00	300.00	32.00
86:00	12:00	300.00	50.00
87:00	01:00	300.00	100.00
94:00	07:00	300.00	100.00
95:00	01:00	300.00	3.00
105:00	10:00	300.00	3.00

Pre-column equilibration:

Volume [µl] : 12.00
Flow [µl / min] : (unspecified)
Max. pressure [Bar] : 800.00

Analytical column equilibration:

Volume [µl] : 3.50
Flow [µl / min] : (unspecified)
Max. pressure [Bar] : 600.00

Autosampler wash:

Flush volume [µl] : 100.00

*Method of Q Exactive HF***Overall method settings****Global Settings**

Use lock masses	best
Lock mass injection	-
Chrom. peak width (FWHM)	20 s

Time

Method duration	105.00 min
-----------------	------------

Customized Tolerances (+/-)

Lock Masses	-
Inclusion	-
Exclusion	-
Neutral Loss	-
Mass Tags	-
Dynamic Exclusion	-

*Experiments***Full MS - SIM****General**

Runtime	0 to 105 min
Polarity	positive
In-source CID	0.0 eV

Full MS - SIM

Microscans	1
Resolution	120,000
AGC target	3e6
Maximum IT	20 ms
Number of scan ranges	1
Scan range	400 to 1000 m/z
Spectrum data type	Profile

DIA**General**

Runtime	0 to 105 min
Polarity	positive
In-source CID	0.0 eV
Default charge state	2

DIA

Microscans	1
Resolution	30,000
AGC target	1e6
Maximum IT	50 ms
Loop count	9
MSX count	1
MSX isochronous ITs	on
Isolation window	23.2 m/z

Isolation offset	0.0 m/z
Fixed first mass	200.0 m/z
(N)CE / stepped (N)CE	nce: 27
Spectrum data type	Profile

Full MS – SIM**General**

Runtime	0 to 105 min
Polarity	positive
In-source CID	0.0 eV

Full MS – SIM

Microscans	1
Resolution	120,000
AGC target	3e6
Maximum IT	20 ms
Number of scan ranges	1
Scan range	400 to 1000 m/z
Spectrum data type	Profile

DIA**General**

Runtime	0 to 105 min
Polarity	positive
In-source CID	0.0 eV
Default charge state	2

DIA

Microscans	1
Resolution	30,000
AGC target	1e6
Maximum IT	50 ms
Loop count	9
MSX count	1
MSX isochronous ITs	on
Isolation window	23.2 m/z
Isolation offset	0.0 m/z
Fixed first mass	200.0 m/z
(N)CE / stepped (N)CE	nce: 27
Spectrum data type	Profile

Full MS – SIM**General**

Runtime	0 to 105 min
Polarity	positive
In-source CID	0.0 eV

Full MS – SIM

Microscans	1
Resolution	120,000
AGC target	3e6
Maximum IT	20 ms
Number of scan ranges	1

Scan range	400 to 1000 m/z
Spectrum data type	Profile

DIA**General**

Runtime	0 to 105 min
Polarity	positive
In-source CID	0.0 eV
Default charge state	2

DIA

Microscans	1
Resolution	30,000
AGC target	1e6
Maximum IT	50 ms
Loop count	9
MSX count	1
MSX isochronous ITs	on
Isolation window	23.2 m/z
Isolation offset	0.0 m/z
Fixed first mass	200.0 m/z
(N)CE / stepped (N)CE	nce: 27
Spectrum data type	Profile

Setup**Tunefiles****General**

Switch Count 0
Base Tunefile C:\Xcalibur\methods\nanosource.mstune

Contact Closure**General**

Used False
Start in Closed True
Switch Count 0

Syringe**General**

Used False
Start in OFF True
Stop at end of run False
Switch Count 0

Pump setup

Syringe type	Hamilton
Flow rate	3.000 µL/min
Inner diameter	2.303 mm
Volume	250 µL

Divert Valve A**General**

Used False
Start in 1-2 True
Switch Count 0

Divert Valve B**General**

Used False
Start in 1-2 True
Switch Count 0

Lock Masses

1 entry

Mass [m/z]	Polarity	Start [min]	End [min]	Comment
445.12003	Positive			

Inclusion List

27 entries

Mass [m/z]	Formula [M]	Species	CS [z]	Polarity	Start [min]	End (N) [min]	CE	MSX ID	Comment
411.00000				Positive					
432.60000				Positive					
454.80000				Positive					
477.00000				Positive					
499.20000				Positive					
521.40000				Positive					
543.60000				Positive					
565.80000				Positive					
588.00000				Positive					
610.20000				Positive					
632.40000				Positive					
654.60000				Positive					
676.80000				Positive					
699.00000				Positive					
721.20000				Positive					
743.40000				Positive					
765.60000				Positive					
787.80000				Positive					
810.00000				Positive					
832.20000				Positive					
854.40000				Positive					
876.60000				Positive					
898.80000				Positive					
921.00000				Positive					
943.20000				Positive					
965.40000				Positive					
987.60000				Positive					

Instrument Method: HCD_90min_grad3_DIA_28iw_TB

Thermo EASY-LC method print for QExactive_HF

Sample pickup:

Volume [µl] : 1.00
Flow [µl / min] : 20.00

Sample loading:

Volume [µl] : 22.00
Flow [µl / min] : (unspecified)
Max. pressure [Bar] : 800.00

Gradient:

Time [mm:ss]	Duration [mm:ss]	Flow [nl/min]	Mixture [%B]
00:00	00:00	300.00	5.00
01:00	01:00	300.00	6.00
51:00	50:00	300.00	27.00
70:00	19:00	300.00	44.00
70:10	00:10	300.00	95.00
80:00	09:50	300.00	95.00
81:00	01:00	300.00	5.00
81:10	00:10	300.00	5.00
90:00	08:50	300.00	5.00

Pre-column equilibration:

Volume [µl] : 12.00
Flow [µl / min] : (unspecified)
Max. pressure [Bar] : 800.00

Analytical column equilibration:

Volume [µl] : 3.50
Flow [µl / min] : (unspecified)
Max. pressure [Bar] : 600.00

Autosampler wash:

Flush volume [µl] : 100.00

*Method of Q Exactive HF***Overall method settings****Global Settings**

Use lock masses	best
Lock mass injection	-
Chrom. peak width (FWHM)	20 s

Time

Method duration	90.00 min
-----------------	-----------

Customized Tolerances (+/-)

Lock Masses	-
Inclusion	-
Exclusion	-
Neutral Loss	-
Mass Tags	-
Dynamic Exclusion	-

*Experiments***Full MS - SIM****General**

Runtime	0 to 90 min
Polarity	positive
In-source CID	0.0 eV

Full MS - SIM

Microscans	1
Resolution	60,000
AGC target	3e6
Maximum IT	40 ms
Number of scan ranges	1
Scan range	400 to 1000 m/z
Spectrum data type	Profile

DIA**General**

Runtime	0 to 90 min
Polarity	positive
In-source CID	0.0 eV
Default charge state	2

DIA

Microscans	1
Resolution	30,000
AGC target	1e6
Maximum IT	40 ms
Loop count	26
MSX count	1
MSX isochronous ITs	on
Isolation window	22.0 m/z

Isolation offset	0.0 m/z
Fixed first mass	200.0 m/z
(N)CE / stepped (N)CE	nce: 27
Spectrum data type	Profile

Setup**Tunefiles****General**

Switch Count 0
Base Tunefile C:\Xcalibur\methods\nanosource.mstune

Contact Closure**General**

Used False
Start in Closed True
Switch Count 0

Syringe**General**

Used False
Start in OFF True
Stop at end of run False
Switch Count 0

Pump setup

Syringe type Hamilton
Flow rate 3.000 µL/min
Inner diameter 2.303 mm
Volume 250 µL

Divert Valve A**General**

Used False
Start in 1-2 True
Switch Count 0

Divert Valve B**General**

Used False
Start in 1-2 True
Switch Count 0

Lock Masses

1 entry

Mass	Polarity	Start	End	Comment
[m/z]		[min]	[min]	
445.12003	Positive			

Inclusion List

28 entries

Mass [m/z]	Formula [M]	Species	CS [z]	Polarity	Start [min]	End (N)CE [min]	MSX ID	Comment
411.00000				Positive				
432.00000				Positive				
453.00000				Positive				
474.00000				Positive				
495.00000				Positive				
516.00000				Positive				
537.00000				Positive				
558.00000				Positive				
579.00000				Positive				
600.00000				Positive				
621.00000				Positive				
642.00000				Positive				
663.00000				Positive				
684.00000				Positive				
705.00000				Positive				
726.00000				Positive				
747.00000				Positive				
768.00000				Positive				
789.00000				Positive				
810.00000				Positive				
831.00000				Positive				
852.00000				Positive				
873.00000				Positive				
894.00000				Positive				
915.00000				Positive				
936.00000				Positive				
957.00000				Positive				
978.00000				Positive				

Instrument Method: HCD_2H5M_T20_HI

Thermo EASY-LC method print for QExactive_HF

Sample pickup:

Volume [µl] : 1.00
Flow [µl / min] : 20.00

Sample loading:

Volume [µl] : 22.00
Flow [µl / min] : (unspecified)
Max. pressure [Bar] : 800.00

Gradient:

Time [mm:ss]	Duration [mm:ss]	Flow [nl/min]	Mixture [%B]
00:00	00:00	300.00	3.00
04:00	04:00	300.00	8.00
06:00	02:00	300.00	10.00
74:00	68:00	300.00	32.00
86:00	12:00	300.00	50.00
87:00	01:00	300.00	100.00
94:00	07:00	300.00	100.00
95:00	01:00	300.00	3.00
105:00	10:00	300.00	3.00

Pre-column equilibration:

Volume [µl] : 12.00
Flow [µl / min] : (unspecified)
Max. pressure [Bar] : 800.00

Analytical column equilibration:

Volume [µl] : 2.00
Flow [µl / min] : (unspecified)
Max. pressure [Bar] : 600.00

Autosampler wash:

Flush volume [µl] : 100.00

*Method of Q Exactive HF***Overall method settings****Global Settings**

Use lock masses	best
Lock mass injection	-
Chrom. peak width (FWHM)	20 s

Time

Method duration	105.00 min
-----------------	------------

Customized Tolerances (+/-)

Lock Masses	-
Inclusion	-
Exclusion	-
Neutral Loss	-
Mass Tags	-
Dynamic Exclusion	-

*Experiment***Full MS / dd-MS² (TopN)****General**

Runtime	0 to 105 min
Polarity	positive
In-source CID	0.0 eV
Default charge state	2
Inclusion	-
Exclusion	-
Tags	-

Full MS

Microscans	1
Resolution	60,000
AGC target	3e6
Maximum IT	32 ms
Number of scan ranges	1
Scan range	350 to 1500 m/z
Spectrum data type	Profile

dd-MS² / dd-SIM

Microscans	1
Resolution	15,000
AGC target	1e5
Maximum IT	50 ms
Loop count	20
MSX count	1
TopN	20
Isolation window	2.0 m/z
Isolation offset	0.0 m/z
Scan range	200 to 2000 m/z
Fixed first mass	110.0 m/z


```

(N)CE / stepped (N)CE                                nce: 26
Spectrum data type                                     Profile
dd Settings
Minimum AGC target                                     1.00e3
Intensity threshold                                     2.0e4
Apex trigger                                           -
Charge exclusion                                       unassigned, 1, 5 - 8, >8
Multiple charge states                                 all
Peptide match                                          preferred
Exclude isotopes                                       on
Dynamic exclusion                                     25.0 s
If idle ..                                           do not pick others

```

Setup**Tunefiles****General**

```

Switch Count    0
Base Tunefile   C:\Xcalibur\methods\nanosource.mstune

```

Contact Closure**General**

```

Used            False
Start in Closed True
Switch Count    0

```

Syringe**General**

```

Used            False
Start in OFF    True
Stop at end of run False
Switch Count    0

```

Pump setup

```

Syringe type    Hamilton
Flow rate       3.000 µL/min
Inner diameter  2.303 mm
Volume          250 µL

```

Divert Valve A**General**

```

Used            False
Start in 1-2    True
Switch Count    0

```

Divert Valve B**General**

```

Used            False
Start in 1-2    True
Switch Count    0

```

1 entry

Mass [m/z]	Polarity	Start [min]	End [min]	Comment
445.12003	Positive			

Instrument Method: HCD_105min_T20_grad3_TB

Thermo EASY-LC method print for QExactive_HF

Sample pickup:

Volume [µl] : 1.00
Flow [µl / min] : 20.00

Sample loading:

Volume [µl] : 22.00
Flow [µl / min] : (unspecified)
Max. pressure [Bar] : 800.00

Gradient:

Time [mm:ss]	Duration [mm:ss]	Flow [nl/min]	Mixture [%B]
00:00	00:00	300.00	4.00
01:00	01:00	300.00	6.00
70:00	69:00	300.00	27.00
85:00	15:00	300.00	44.00
85:10	00:10	300.00	95.00
95:00	09:50	300.00	95.00
95:10	00:10	300.00	5.00
104:50	09:40	300.00	5.00
105:00	00:10	300.00	5.00

Pre-column equilibration:

Volume [µl] : 12.00
Flow [µl / min] : (unspecified)
Max. pressure [Bar] : 800.00

Analytical column equilibration:

Volume [µl] : 2.00
Flow [µl / min] : (unspecified)
Max. pressure [Bar] : 600.00

Autosampler wash:

Flush volume [µl] : 100.00

*Method of Q Exactive HF***Overall method settings****Global Settings**

Use lock masses	best
Lock mass injection	-
Chrom. peak width (FWHM)	20 s

Time

Method duration	105.00 min
-----------------	------------

Customized Tolerances (+/-)

Lock Masses	-
Inclusion	-
Exclusion	-
Neutral Loss	-
Mass Tags	-
Dynamic Exclusion	-

*Experiment***Full MS / dd-MS² (TopN)****General**

Runtime	0 to 105 min
Polarity	positive
In-source CID	0.0 eV
Default charge state	2
Inclusion	-
Exclusion	-
Tags	-

Full MS

Microscans	1
Resolution	60,000
AGC target	3e6
Maximum IT	32 ms
Number of scan ranges	1
Scan range	350 to 1500 m/z
Spectrum data type	Profile

dd-MS² / dd-SIM

Microscans	1
Resolution	15,000
AGC target	1e5
Maximum IT	50 ms
Loop count	20
MSX count	1
TopN	20
Isolation window	2.0 m/z
Isolation offset	0.0 m/z
Scan range	200 to 2000 m/z
Fixed first mass	110.0 m/z

```

(N)CE / stepped (N)CE                                nce: 26
Spectrum data type                                     Profile
dd Settings
Minimum AGC target                                     1.00e3
Intensity threshold                                     2.0e4
Apex trigger                                           -
Charge exclusion                                       unassigned, 1, 5 - 8, >8
Multiple charge states                                 all
Peptide match                                          preferred
Exclude isotopes                                       on
Dynamic exclusion                                     25.0 s
If idle ..                                           do not pick others

```

Setup**Tunefiles****General**

```

Switch Count    0
Base Tunefile   C:\Xcalibur\methods\nanosource.mstune

```

Contact Closure**General**

```

Used            False
Start in Closed True
Switch Count    0

```

Syringe**General**

```

Used            False
Start in OFF    True
Stop at end of run False
Switch Count    0

```

Pump setup

```

Syringe type    Hamilton
Flow rate       3.000 µL/min
Inner diameter  2.303 mm
Volume          250 µL

```

Divert Valve A**General**

```

Used            False
Start in 1-2    True
Switch Count    0

```

Divert Valve B**General**

```

Used            False
Start in 1-2    True
Switch Count    0

```

1 entry

Mass [m/z]	Polarity	Start [min]	End [min]	Comment
445.12003	Positive			